

Whole-Genus Sequencing: 300 Aspergilli

Alan Kuo^{*1}, Robert Riley¹, Alicia Clum¹, Asaf Salamov¹, Scott Baker², Blake Simmons²,
Mikael Andersen³, Igor Grigoriev¹

¹DOE Joint Genome Inst., Walnut Creek, CA,

²Joint BioEnergy Inst., Emeryville, CA

³Technical U. of Denmark, Kgs. Lyngby, Denmark

March 2015

The work conducted by the U.S. Department of Energy Joint Genome Institute is supported
by the Office of Science of the U.S. Department of Energy under Contract No. DE-AC02-
05CH11231

LBNL-178198

DISCLAIMER

This document was prepared as an account of work sponsored by the United States Government. While this document is believed to contain correct information, neither the United States Government nor any agency thereof, nor The Regents of the University of California, nor any of their employees, makes any warranty, express or implied, or assumes any legal responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by its trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof, or The Regents of the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof or The Regents of the University of California.


ABSTRACT

Aspergillus is a ubiquitous and phenotypically diverse genus of filamentous Ascomycota, many of which play key roles as fermenters in food production, platforms for biotechnology and industrial production of enzymes and chemicals, plant and opportunistic animal pathogens, and agents of agricultural toxigenesis and biomass conversion for bioenergy. As part of a DOE Joint BioEnergy Institute initiative to characterize the entire genus, the JGI plans to sequence, assemble, and annotate the genomes of each of the ~300 species of the genus *Aspergillus*. To accomplish this massive task in a timely manner without sacrificing quality, we have sought to streamline our existing processes as well as explore alternative technologies, especially assembly and annotation of long PacBio sequencing reads. Over the past year we have released on MycoCosm the genomes of an additional 24 *Aspergillus* sp. with preliminary analyses of their phylogenies, secretomes, and secondary metabolism. The next tranche of 92 species is expected soon.

Why whole-genus sequencing?

The 300+ species of Aspergilli have been identified and classified based on morphology, and on profiles of secondary metabolites, of which *Aspergillus* sp. produce a greater variety than any other fungi.

Fig. 1: Two of the 24 species used in this study, each with distinct morphology and secondary metabolism on 6 standard growth media.




Advances in genomics allow increasing numbers of aspergilli to be analyzed, and now we are in the process of sequencing every species in the genus. An entirely-sequenced genus will provide a comprehensive platform for taxonomy and phylogeny, as well as functional studies including but not limited to secreted proteins and secondary metabolites.

To optimize sequencing strategies for such a grand scale, JGI is exploring ways to streamline our standard processes. We are testing on some of the first 24 *Aspergillus* sp. we have received for this project. We can also provide some preliminary results for these 24 genomes.

Credits


Ellen Lyhne (DTU)
cultures, photos
Martin Kogle (DTU)
DNAs, RNAs
Matt Nolan (JGI)
genome assembly
Anna Lipzen (JGI)
transcriptome assembly
Frank Korzeniewski (JGI)
SMC pipeline
Roman Nikitin (JGI)
SMC portal



Explore alternative sequencing technology

For 6 *Aspergillus* sp. we did both PacBio sequencing + assembly, as well as JGI standard Illumina sequencing + assembly, followed by the standard JGI Annotation Pipeline. We then compared the resulting assemblies and annotations.


Fig. 2: Representative results with *A. eucalypticola*. PacBio did not provide a significant advantage over Illumina,



Explore reductions in RNASeq

For 4 *Aspergillus* sp. we reduced the number of RNASeq reads used for JGI standard transcriptomic assembly and for JGI Annotation Pipeline. We then compared the resulting annotations.

Fig. 3: Representative results with *A. steynii*. Quality began to fall at 1/6th the standard no. of reads.



Preliminary results

Fig. 4A: ML species tree using 2k conserved genes. Most subgenera are confirmed.


Fig. 4B: Predicted gene complements and secreted proteins.

Species	Genome (Mbp)	1000s genes	# secreted proteins
Aspjap1 ↔ <i>A. japonicus</i>	36.1	12.0	1374
Aspvio1 ↔ <i>A. violaceofuscus</i>	36.0	12.1	1377
Aspind1 ↔ <i>A. indologenus</i>	38.6	12.1	1382
Aspuva1 ↔ <i>A. uvarum</i>	35.9	12.0	1367
Aspbru1 ↔ <i>A. brunneoviolaceus</i>	37.5	12.1	1396
Aspfij1 ↔ <i>A. fijiensis</i>	36.5	12.0	1366
Aspacu1 ↔ <i>A. aculeatinus</i>	36.5	12.0	1387
Aspsac1 ↔ <i>A. saccharolyticus</i>	31.1	10.1	1061
Asphom1 ↔ <i>A. homomorphus</i>	34.1	11.4	1205
Aspeuc1 ↔ <i>A. eucalypticola</i>	34.8	11.9	1310
Asppip1 ↔ <i>A. piperis</i>	35.3	12.1	1373
Aspneo1 ↔ <i>A. neoniger</i>	35.4	11.9	1318
Aspcos1 ↔ <i>A. costaricensis</i>	36.9	12.0	1330
Aspvad1 ↔ <i>A. vadensis</i>	35.7	12.1	1342
Aspni7 ↔ <i>A. niger</i>	34.9	11.9	1300
Asplac1 ↔ <i>A. lacticoffeatus</i>	35.9	13.1	1388
Aspibe1 ↔ <i>A. ibericus</i>	33.4	11.7	1321
Aspscle1 ↔ <i>A. sclerotiiicarbonarius</i>	37.6	12.6	1401
Aspscl1 ↔ <i>A. sclerotioniger</i>	36.7	12.3	1295
Aspell1 ↔ <i>A. ellipticus</i>	42.9	12.9	1351
Asphet1 ↔ <i>A. heteromorphus</i>	35.6	11.1	1183
Aspnd1 ↔ <i>A. nidulans</i>	30.5	10.7	1142
Aspoch1 ↔ <i>A. ochraceoroseus</i>	27.7	8.9	839
Aspcam1 ↔ <i>A. campestris</i>	28.3	9.8	1058
Aspste1 ↔ <i>A. steynii</i>	37.8	13.2	1400
Aspor1 ↔ <i>A. oryzae</i>	37.9	12.0	1333
Aspfu1 ↔ <i>A. fumigatus</i>	29.4	9.8	1026
Aspnov1 ↔ <i>A. novofumigatus</i>	32.4	11.5	1299
Pench1 ↔ <i>P. chrysogenum</i>	31.3	11.4	1083

Fig. 4C: Predicted secondary metabolite clusters (SMCs). Statistically significant gains are blue, losses are red.

	PKS	PKS-Like	NRPS	NRPS-Like	HYBRID	DMAT	TC
19	3	21	22	1	4	1	
17	4	19	20	3	6	3	
27	3	17	24	5	6	2	
21	3	17	21	3	4	2	
24	3	21	21	2	4	4	
26	3	21	25	3	5	4	
28	3	19	20	2	6	3	
11	6	14	15	2	3	0	
24	5	16	22	3	4	4	
26	4	11	19	6	2	6	
30	5	14	20	6	2	8	
29	6	13	16	4	3	9	
32	7	17	21	5	3	7	
27	5	13	22	5	2	4	
27	4	16	16	9	2	6	
23	6	16	16	9	2	7	
15	6	7	18	6	2	3	
25	6	14	22	6	3	3	
21	5	9	21	5	3	5	
28	7	15	14	2	2	4	
16	3	12	16	3	4	3	
22	5	9	13	1	5	1	
10	2	6	3	2	4	1	
12	2	11	11	5	5	2	
29	3	17	22	9	9	4	
23	5	14	15	2	9	2	
14	2	9	7	0	2	0	
24	2	16	5	3	4	2	
18	3	9	14	2	0	4	

Fig. 4D: Some SMCs are conserved across the genus



Conclusions, so far

- Methodological optimization, so far**
- PacBio sequencing of *Aspergillus* genomes was not substantially better than Illumina.
 - RNASeq was reduced 6-fold without adverse effects on *Aspergillus* annotation.
 - SMC pipeline and portal were developed and tested.
- Biological results, so far**
- 24 *Aspergillus* sp. across multiple subgenera were assembled and annotated.
 - Tentative phylogeny was inferred using whole-genome methods.
 - Secreted proteins were identified.
 - SMCs were cataloged and clade-specific gains and losses were identified.

What's next

- More streamlining of methods**
- One example: test pipelines with fewer gene predictors.
- More analyses**
- Target other functionally significant proteins: proteases, transporters, CAZymes.
 - Identify gains and losses in the above categories of genes
- More genomes**
- Next tranche of 92 species is in sequencing stage.